

%	n_D^{20}	%	n_D^{20}	%	n_D^{20}	%	n_D^{20}	%	n_D^{20}
0	1.33299	20	1.36384	40	1.39986	60	1.44193	80	1.49071
1	1.33442	21	1.36551	41	1.40181	61	1.44420	81	1.49333
2	1.33586	22	1.36720	42	1.40378	62	1.44650	82	1.49597
3	1.33732	23	1.36889	43	1.40576	63	1.44881	83	1.49862
4	1.33879	24	1.37060	44	1.40776	64	1.45113	84	1.50129
5	1.34026	25	1.37233	45	1.40978	65	1.45348	85	1.50398
6	1.34175	26	1.37406	46	1.41181	66	1.45584	86	1.5067
7	1.34325	27	1.37582	47	1.41385	67	1.45822	87	1.5094
8	1.34477	28	1.37758	48	1.41592	68	1.46061	88	1.5122
9	1.34629	29	1.37936	49	1.41799	69	1.46303	89	1.5149
10	1.34782	30	1.38115	50	1.42009	70	1.46546	90	1.5177
11	1.34937	31	1.38296	51	1.42220	71	1.46790	91	1.5205
12	1.35093	32	1.38478	52	1.42432	72	1.47037	92	1.5234
13	1.35250	33	1.38661	53	1.42647	73	1.47285	93	1.5262
14	1.35408	34	1.38846	54	1.42863	74	1.47535	94	1.5291
15	1.35568	35	1.39032	55	1.43080	75	1.47787	95	1.5320
16	1.35729	36	1.39220	56	1.43299	76	1.48040		
17	1.35891	37	1.39409	57	1.43520	77	1.48295		
18	1.36054	38	1.39600	58	1.43743	78	1.48552		
19	1.36218	39	1.39792	59	1.43967	79	1.48811		